

**Application**

Ergonomically designed flexible belt to prevent it slipping down. Detachable Velcro and separate belt for various applications.

**Combined dual heads**



**Separated head**



**Design**



- ▶ Alternative colorful design
- ▶ Smaller and compact size
- ▶ Powerful movement

**BUHEUNG MEDICAL Co., Ltd.**  
 336-30 Doksan1-dong Geumcheon-gu Seoul 153-814 Republic of Korea  
 www.buheungmedical.com  
 overseas@buheung.co.kr  
 Tel: +82 2 838 5333 Fax: +82 2 838 5334

Multi-Functional Exercise Machine  
 "FreeM plus"



FreeM+

All in one!

- ✓ Asymmetrical Dual Vibration
- ✓ Optimized Combined EMS
- ✓ Body Fat Check
- ✓ Facial Lifting Function
- ✓ Music Therapy

**BUHEUNG MEDICAL Co., Ltd.**

**MK-307 Multi-Functional Exercise Machine**

FreeM+

**Optimized multiple effects with vibration, EMS, facial lifting and music therapy.**


- Asymmetrical movement with Dual Heads safely transfers the vibration waves deeply to the body.
- Optimized combined frequencies (ICT / Middle Frequency / Low Frequency / Electric Acupuncture) burn the fat and glycogen in the muscles due to muscle contraction.
- Facial lifting function using wireless remote control.
- Music can affect the body in many health-promoting ways. This is the basis for a growing field of therapy known as music therapy.

**Multiple Effects with Vibration and EMS**



**EMS + VIBRATION**  
 TAR Motors (Two-way Asymmetrical Rotating Motors) and Frequency Technology maximize its slimming effect whilst asymmetrical CYCLONE movement with middle and low frequency stimulation is simultaneously penetrating into the diet points.

**VIBRATION**




**Asymmetry movement with TAR motor**  
 TAR Motor System (Two-way Asymmetry Rotating Motors) provides you dynamic and effective CYCLONE massage for your slimming and health care. This system accelerates the blood and lymph flow helping to reduce fatty deposits and cellulite.

**EMS (ELECTRIC MUSCLE STIMULATION)**


**MOST EFFECTIVE FREQUENCIES**  
 Reduction in patient tolerance effect using Frequency Technology with various kinds of wave characteristics. A chemical reaction results from muscle contraction and it burns fat and glycogen in the muscles, which enables you to improve your body contours. With frequencies 1~30Hz, stimulation can have an acupuncture effect on the troubled parts of your body.

**DEEP FREQUENCY PENETRATION**



Frequencies consisting mostly of middle frequencies have improved penetrating power deep into your body.

**MORE EFFECTIVE ICT (Interferential Current Therapy)**




Using ICT, which is generated by two different kinds of middle frequencies, FreeM+ reduces pain related to muscles, the skeleton and the nervous system.

**BODY FAT MEASUREMENT**

▶ With the body fat sensor on the Rolling contact pads you can measure the body fat value in any area of the body before and after application of FreeM+.

The user easily identifies the application area before exercising and then checks the exercise result with the body fat measurement.

The body fat value consists of numerals from 001 to 999.



**MUSIC THERAPY**

Music physically changes blood pressure, pulse rate, respiration, skin reaction, brain wave patterns and muscle reactions and psychologically achieves stress relief benefits.



Optimized music therapy multiplies the slimming effect with Vibration and EMS functions.

**Facial Lifting**

Facial skin care with remote control whilst using the belt massager on any area of your body provides additional time to care for other areas of your body. With middle frequencies, stimulation can help a facial lifting effect.

- ◊ Wireless operation with remote control.
- ◊ Easy to choose the treatment area with the remote control.



**TOUCH PANEL**



- ◊ Touch Quick control panel on the main head employs advanced technology.

**Rolling electrodes**



- ▶ With the help of 8 rolling electrodes for blood circulation, and skin elasticity, it is easy to move the FreeM+ heads.
- ▶ The electrode rolls, stimulating with various kinds of frequencies.
- ▶ The distance between the Rolling Electrode pads is specially designed to stimulate acupuncture points.
- ▶ Nano-Silver Technology for Antimicrobial Efficacy

**SPECIFICATIONS**

- ◆ Model : FreeM+ (MK-307)
- ◆ Input : 100 ~ 240V AC / 50Hz, 60Hz
- ◆ Output : 24V DC1.9A

