

# SHELLA™

## MK-902



SHELLA



*ANTI-Aging  
Collagen Activation  
Blood circulation  
Cosmetics Penetration  
Whitening  
Cleansing  
iontophoresis*



### Features of SHELLA™ MK-902

- Heating & Cooling at the same time, different place
- Electro Collagen Stimulation
- (-) iontophoresis
- (+) iontophoresis
- Ergonomic & Beautiful Design
- Portable Skin care device



## SHELLA MK-902

20℃ electro frequency skin science

Model No. / Name	MK-902 SHELLA
Function	Heating Cooling Electro Medium Frequency Option) (-) iontophoresis (+) iontophoresis
Material	Titanium, PC, Silicon, etc.
Power	3.7V / 1500mA mini USB rechargeable battery
Manufacturer / Origin	Buheung Medical / Made in Korea
Size / Weight	123 x 76 x 29 (mm) / 135g

#### ► Effect ◀

- Collagen Expression increase
- Cosmetics Penetration increase
- Blood Circulation
- Cleansing
- Bacterial Skin trouble improvement
- Skin condition improvement
- Edema improvement

#### ► Maximize Skin Condition Improvement! ◀

- Heating & Cooling at the same time
- Micro current Mid -frequency for Skin
  - Collagen Activation
  - No discomfort
- (-) & (+) iontophoresis
  - (-) iontophoresis supplies cosmetics into skin
  - (+) iontophoresis pull dusts from skin. (Cleansing)
- Simple & Easy Usage
- Ergonomic & Beautiful Design
  - For All part of face.
  - Comfortable Grip, portable size, Beautiful Design.
  - Much Less time, Maximized Effect!
- Safety
  - 3 processes temperature control system.
  - Titanium Plate
  - Automatic Power off system

#### ► Effect & Feature for each Function ◀

- Heating & Cooling
  - Heating & Cooling At the same time  
Front Plate: Heating / Back Plate: Cooling
  - Effective & safe temperature control  
Heating: 35~41°C / Cooling: 15~20°C
  - Effects
    - Blood Circulation improvement
    - Collagen Expression increase
    - Muscle, Pores stimulation
    - Bacterial Skin trouble improvement
    - Edema improvement
- Electro Skin Therapy Current
  - Effect with Heating & Cooling
  - No discomfort Feeling
  - Effect
    - Cell & Collagen Activation
    - Skin condition improvement
    - Front Plate Current  
Cosmetics penetration increase  
Collagen Expression increase
    - Back Plate Current  
Lifting Effect
- (-) iontophoresis
  - Cosmetics Penetration increase
- (+) iontophoresis
  - Cleansing

#### ► Effective Usage ◀

- Heating, Cooling & Micro current Mid - Frequency
  - ① 10 sec. Heating massage and 5 sec. Cooling Massage on your specific face part repeatedly. (2:1 Method)
  - ② With creams, spread cream with heating and cooling plate as ① Method.
  - ③ With Maskpack, Heating and cooling massage on the maskpack as ① method.
  - ④ When Edema, 10 Sec. Heating and 10 Sec. Cooling massage repeatedly.
- (-) iontophoresis
  - ① Use the cream, skin, lotion or serums on the facial skin.
  - ② Spread cosmetics with front metal plate.
- (+) iontophoresis
  - ① Use the cleansing cream on the facial skin.
  - ② Spread cleansing cream with front metal plate.